## INTENDED USE

DOAC-Stop can be used to efficiently extract all types of Direct Oral Anti-Coagulants (DOACs) including dabigatran, edoxaban, betrixaban, rivaroxaban, apixaban and also argatroban from test plasmas with minimal effect on plasma proteins involved in the clotting mechanism. This product is intended for use as an accessory IVD, and it is CE marked.

## INTRODUCTION

Therapeutic uses of DOACs are increasing. DOACs are known to interfere with almost all clotting tests to varying degrees and sometimes patients who need to be tested for underlying coagulation defects may also be on DOACs.

Specific antidotes for individual DOACs are being developed for therapeutic use but are not yet widely available for laboratory use. DOAC-Stop is the first general agent available for solving the diagnostic problems associated with DOACs. After treatment with DOAC-Stop, plasma samples can be analysed for underlying coagulation defects such as factor deficiencies, heparin, lupus anticoagulant or other interfering antibodies [1].

## CONTENTS OF PRODUCT

| Product Code | Pack Size   |
|--------------|-------------|
| HX9904-100   | 100 tablets |
| HX9904-50    | 50 tablets  |

## LIMITATIONS

DOAC-Stop may extract low molecular weight compounds resembling DOACs from plasmas. Those affecting coagulation (but administered intravenously, not orally) include argatroban, aprotinin, protamine, hirudin analogues, polymyxin and related cationic compounds [2].

## PRECAUTIONS

DOAC Stop is intended for use with plasma samples suspected to contain DOACs. If test results are unchanged by DOAC Stop and DOACs are still suspected to be present, apply appropriate chromogenic anti factor Xa or anti thrombin assays to obtain specific DOAC results. Contact your distributor or manufacturer for technical support.

Store at 2-8°C. Keep Dry. Do not use after the expiry date indicated on the label. Treat all clinical material as potentially infectious and dispose in accordance with local operating regulations. For further information, please refer to Safety Data Sheet and Product Information.

## INSTRUCTIONS FOR USE

**Sample preparation:** DOAC-Stop has been developed for use with citrated plasmas. Follow your validated laboratory procedures for preparing test plasma. Apply the citrated plasma in the procedure below to remove any DOAC if present.

Method for DOAC removal from a test plasma:

|    |                                                                                                                                                            |
|----|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1. | Add 1.0mL of citrated test plasma to be treated to a plastic centrifuge tube. <i>Ideally 1.0mL plasma but a range from 0.5 ml to 1.5 ml is acceptable.</i> |
| 2. | Add 1 DOAC-Stop mini-tab and mix gently until it has dispersed. Mix for a further 5-10 minutes.                                                            |
| 3. | Centrifuge down the particulate matter (for example, 5 minutes at 2000g or 1 minute at 7000g in a microfuge). Do not stop centrifuging too quickly.        |
| 4. | The supernatant plasma, now depleted of DOACs can be used directly for testing or can be transferred into a separate plastic tube for frozen storage.      |

## APPLICATION

Plasmas treated with DOAC-Stop can be used for reliable routine assessment of known blood coagulation parameters [3-7]. The degree of shortening induced by DOAC-Stop in a clotting time test depends on the concentration and type of DOAC as well as the on the underlying plasma itself. Clotting time results with DOAC sensitive tests such as dRVVT and APTT tests, may be expressed as ratios of test result before (B) to that obtained after (A) DOAC Stop treatment. The B/A result ranges upward from 1.0 and correlates positively with DOAC concentration depending on which DOAC and test is used [8].

## PERFORMANCE CHARACTERISTICS

One mini-tab of DOAC Stop in 1.0mL of normal plasma spiked with 500ng/ml of dabigatran, edoxaban, betrixaban, rivaroxaban or apixaban will remove more than 95% of the DOAC within 5 minutes. There was no effect on the baseline APTTs after 3 hours further incubation. In 92 test plasma samples covering a range of abnormalities, 89 gave definitive results and 3 were equivocal after DOAC-Stop in view of patient complexity [1].

A recent study on 135 DOAC-treated plasmas showed that DOAC Stop was mostly effective in overcoming the effects of dabigatran in 40 cases, apixaban in 38 cases, rivaroxaban 42 cases and edoxaban 15 cases on several thrombophilia screening tests. The incidence of false positive lupus anticoagulants fell from approximately 30% to zero [3].

## INDEMNITY NOTICE

DOAC-Stop is intended to be used with plasma samples containing DOACs. Follow procedures and refer to precautions that may affect the stated or implied claims and performance of this product. Cellabs, Haematex and its agents or distributors are not liable for damages.

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## REFERENCES

- [1] Exner T, et al. "Simple method for removing DOACs from plasma samples" Throm Res. 2018; 16: 1028-39.
- [2] Exner T, et al. "Effect of an activated charcoal product (DOAC Stop™) intended for extracting DOACs on various other APTT-prolonging agents". Clin Chem Lab Med. 2019; 57: 690-696.
- [3] Favresse J, et al. "Evaluation of the DOAC Stop procedure to overcome the effect of DOACs on several thrombophilia screening tests". TH Open, 2018; 2: e202-e209.
- [4] Platton S, Hunt C. "Influence of DOAC Stop on coagulation assays in samples from patients on rivaroxaban or apixaban". Int J Lab Haematol. 2019; 41: 227-233.
- [5] Jacquemin M, et al. "The adsorption of dabigatran is as efficient as addition of Idarucizumab to neutralize the drug in routine coagulation assays. Int J Lab Hematol. 2018; 40: 442-447".
- [6] Favalaro EJ, et al. "Neutralising rivaroxaban-induced interference in laboratory testing for lupus anticoagulants (LA): A comparative study using DOAC Stop and Andaxanet alpha". Throm Res. 2019; 180:10-19.
- [7] Kopatz WF, et al. "Use of DOAC Stop for the elimination of anticoagulants in the thrombin generation assay". Throm Res. 2018; 170: 97-101.
- [8] Exner T, et al "Clotting tests correlate better with DOAC concentrations when expressed as a "Correction Ratio"; results before/after extraction with the DOAC Stop reagent. Throm Res.2019; 179: 69-72.
- [9] Exner T, et al "The effect of DOACs on laboratory tests and their removal to limit interference in functional assays. Int J Lab Hematol. 2020; 42 (Suppl.1): 41-48.

## EXPLANATION OF SYMBOLS EXPLICATION DES SYMBOLES ERKLAERUNG DER SYMBOLE SPIEGAZIONE DEI SIMBOLI

|                                                                                     |                                                                                                                                                      |                                                                                     |                                                                                                                                                                                 |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Consult instructions for use<br>Consultez la notice du produit<br>Gebrauchsanleitung konsultieren<br>Consultare le istruzioni per l'uso              |  | Catalogue number<br>Numéro de catalogue<br>Katalognummer<br>Codice del Prodotto                                                                                                 |
|  | <i>In Vitro</i> Diagnostic Medical Device<br>Produit de diagnostic médical in-vitro<br><i>In Vitro</i> -Diagnostikum<br>Per uso Diagnostico In Vitro |  | Contains sufficient for <n> tests<br>Contient suffisamment de réactifs pour <x> tests<br>Enthält genügend für <n> Tests<br>Contiene la quantità sufficiente per <n> test        |
|  | Temperature limitation<br>Limites de température<br>Temperaturbegrenzung<br>Limitazioni di Temperatura                                               |  | Manufacturer:                      Hersteller:<br>Fabricant :                      Produttore:<br><br>Cellabs, Unit 7/27 Dale St. Brookvale, NSW 2100, Australia                |
|  | Keep Dry<br>Conserver au sec<br>Trocken halten<br>Mantenere in luogo asciutto                                                                        |  | EC Representative                      EG-Vertreter<br>Représentant CE                      Legale Rappresentante<br><br>WMDE B. V, Bergerweg 18, 6085 AT Horn, The Netherlands |
|  | Batch code<br>Lot de produit<br>Batch code<br>Numero di Lotto                                                                                        | Insert<br>version                                                                   | <br><br>LH1.11 28 Sept 2021                                                                  |
|  | Expiry Date<br>Date d'expiration<br>Verfallsdatum<br>Data di Scadenza                                                                                |  |                                                                                                                                                                                 |